## The Delphion Integrated View

 Buy N w: PDF | More choices...

 Vi w: Expand Details | INPADOC | Jump to: Top
 ✓ Go to: Derwent
 ✓ Email this to a friend

**PTitle:** WO0068692A1: A METHOD OF DETECTING AN ANALYTE USING

SEMICONDUCTOR NANOCRYSTALS[French]

P Derwent Title: Detection of analytes using semi-conductor nanocrystals which are more

robust than organic fluorescent dyes and which can be made to have

characteristic spectral emissions [Derwent Record]

**VO World Intellectual Property Organization (WIPO)** 

**%**Kind: A1 Publ.of the Int.Appl. with Int.search report i

PInventor: BRUCHEZ, Marcel, P.; 32392 Derby Street, Union City, CA 94587,

United States of America

DANIELS, R., Hugh; 136 Seale Avenue, Palo Alto, CA 94301, United

Kingdom

EMPEDOCLES, Stephen, A.; 2507 Mardell Way, Mountain View, CA

94043, United States of America

PHILLIPS, Vince, A.; 863 Lewis Avenue, Sunnyvale, CA 94086, United

States of America

WONG, Edith, Y.; 22 Volteraa Court, Danville, CA 94526, United States

of America

ZEHNDER, Donald, A.; 1398 Orange Avenue, San Carlos, CA 94070,

United States of America

SAssignee: QUANTUM DOT CORPORATION, 4030 Fabian Way, Palo Alto, CA

94303, United States of America

News, Profiles, Stocks and More about this company

Published / Filed: 2000-11-16 / 2000-05-05

PApplication WO2000US0012227

Number:

**@IPC Code: G01N 33/58**; G01N 33/533; G01N 33/542;

**PECLA Code: G01N33/533**; G01N33/543D; G01N33/58J;

Priority Number: 1999-05-07 **US1999000133084P** 

PAbstract:

The use of semiconductor nanocrystals as detectable labels in various chemical and biological applications is disclosed. The methods find use for detecting a single analyte, as well as multiple analytes by using more than one semiconductor nanocrystal as a detectable label, each of which emits at a distinct wavelength.

[French]

Attorney, Agent ROBINS, Roberta, L.;

or Firm:

**VINPADOC** Show legal status actions


I status actions

Buy Now: Family Legal Status Report

Related Applications:

Legal Status:

Applicati n NumberFiledPatentPub. DateTitleUS1999000133084P1999-05-07



P Designated

Country:

AE AG AL AM AT AU AZ BA BB BG BR BY CA CH CN CR CU CZ DE DK DM DZ EE ES FI GB GD GE GH GM HR HU ID IL IN IS JP KE KG KP KR KZ LC LK LR LS LT LU LV MA MD MG MK MN MW

MX NO NZ PL PT RO RU SD SE SG SI SK SL TJ TM TR TT TZ UA UG US UZ VN YU ZA ZW, Eur pean patent: AT BE CH CY DE DK ES FI FR GB GR IE IT LU MC NL PT SE, OAPI pat nt: BF BJ CF CG CI CM GA GN GW ML MR NE SN TD TG, ARIPO patent: GH GM KE LS MW SD SL SZ

TZ UG ZW, Eurasian patent: AM AZ BY KG KZ MD RU TJ TM

Family: Show 9 known family members

PDescription: Expand full description

**± A METHOD OF DETECTING AN ANALYTE USING** 

SEMICONDUCTOR NANOCRYSTALS

**+ FIELD OF THE INVENTION** 

**±** BRIEF DESCRIPTION OF THE DRAWINGS

**+ DETAILED DESCRIPTION OF THE INVENTION** 

First Claim: Show all claims

1.A method of detecting one or more target analy-tes in a sample containing or suspected of containing the one or more analytes, comprising the steps of(a) providing the sample on a solid support; (b) combining said sample with a semiconductor nanocrystal conjugate, wherein said combining is performed under conditions that allow formation of a complex comprising said conjugate and said analyte, when present; (c) removing any unbound conjugate; (d) detecting the presence of the complex, if present, by monitoring a spectral emission mediated by the semiconductor nanocrystal in the complex, wherein the emission indicates the presence of one or

more target analytes in the sample. †

©Other Abstract Info:

CHEMABS 133(25)346793J CHEMABS 133(25)346793J DERABS C2001-061110 DERABS C2001-

061110



Powered by







Nominate this for the Gallery...



© 1997-2003 Thomson Delphion

Research Subscriptions | Privacy Policy | Terms & Conditions | Site Map | Contact Us | Help